# nedis

PESC110DCT

MANUAL (p. 2)

ANLEITUNG (S. 4)

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MANUALE (p. 10)

MANUAL DE USO (p. 12)

MANUAL (p. 14)

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KÄYTTÖOHJE (s. 18)

BRUKSANVISNING (s. 20)

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MANUAL DE UTILIZARE (p. 24)

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## Personal Scale

## **Advice for use and care:**

Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surfaces. Always weigh yourself with the same scale placed on the same surface. Uneven floors may affect the reading.

The surface of the scale will be slippery if wet. Keep it dry!

Weigh yourself without clothing and footwear, before meals and always at the same time of day. Keep the scale in a dry place to protect the electronic components. Clean after use with a lightly damp cloth, do not use solvents or immerse the product in water. Avoid overloading the scale, otherwise, permanent damage may occur. A user who has any metal material inserted into his/her body may get differences in measurement of body fat reading. Pacemaker users are advised against using this device.

The scale is designed for household use only, not to be used for commercial purposes.

## **Operating instructions:**

1. Put the scale on a hard flat surface (avoiding carpets or any soft surfaces).
2. Switch on the scale by a slight tap on the glass platform near the LCD display. Then wait a few seconds, the following will appear: 0.0 kg.

**Note:** The scale will not work properly if you step on the scale before it shows 0.0. If this happens please start from the beginning of step 2.

3. Gently step on the scale. Parallel your feet on the scale without moving and wait whilst the scale computer weighs.

When you step off the scale, it will automatically switch off. (If no weighing is carried out while the zero is being displayed, the scale will automatically switch themselves off after approximately 6 seconds).

When the scale is overloaded (over 180 kg/396 lb), the following will appear on the screen: Err.

When the battery needs replacement, the following will appear on the screen: LO.

## **Specifications:**

- Equipped with a high precision "strain-gauge" sensor.
- Capacity: 180 kg/396 lb
- Graduation: 0.1 kg/0.2 lb/0.25 lb
- LCD display
- Foot-tap switch
- Automatic zero resetting
- Automatic switch off
- Overload indicator
- Low battery indicator

Lithium battery is included (CR2032). Please remove the insulation sheet from the battery contacts before using.

**Safety precautions:**

This product should ONLY be opened by an authorized technician when service is required. Do not expose the product to water or moisture..

**Maintenance:**

Clean only with a dry cloth.

Do not use cleaning solvents or abrasives.

**Warranty:**

No guarantee or liability can be accepted for any changes and modifications of the product or damage caused due to incorrect use of this product.

**General:**

Designs and specifications are subject to change without notice. All logos, brands and product names are trademarks or registered trademarks of their respective holders and are hereby recognized as such.

**Attention:**

This product is marked with this symbol. It means that used electrical and electronic products should not be mixed with general household waste. There is a separate collections system for these products.

This product has been manufactured and supplied in compliance with all relevant regulations and directives, valid for all member states of the European Union. It also complies to all applicable specifications and regulations in the country of sales.

Formal documentation is available upon request. This includes, but is not limited to: Declaration of Conformity (and product identity), Material Safety Data Sheet, product test report.

Please contact our customer service desk for support:

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NEDIS B.V., De Tweeling 28, 5215 MC 's-Hertogenbosch, THE NETHERLANDS (HOLLANDA)