

Top-quality, mercury free hearing aid batteries, specially designed for the latest high-performance devices. Optimized for wireless, Bluetooth and MFI devices. Active Core Technology: Rayovac Active Core Technology provides a high quality battery by boosting reliability over the battery's lifespan and reducing hearing aid intermittency issues.



Information

- Guaranteed storage time of 4 years
- 110 Years of experience
- Clear and visible size designation
- Icons to indicate usage and functions
- Strong color coding to guide consumers
- Full product guarantee
- Convenient, easy to handle packaging, clear color coding

Article Details

Ordernumber	04610 945 416
Customs Tariff Number	85066000
Country Of Origin	United Kingdom
IEC	PR70
Others	V10
System	Primary Zinc Air
Voltage	1.45
Consumer Unit	6
Packaging Information	Single Blister

GTIN / EAN



Logistical Data

	Content (SU)	GTIN / EAN	Facing (mm)	Depth (mm)	Height (mm)	Weight (kg)
VBE (Selling unit)	1	50 00252 003212	50.0	9.0	108.0	0.006
SK2 (Inner carton 2)	10	50 00252 003250	110.0	60.0	53.0	0.075
SK1 (Inner carton 1)	100	50 00252 003298	120.0	315.0	120.0	0.828
VSK (Master carton)	1000	50 00252 003335	330.0	615.0	260.0	8.700
PAL (Pallet)	16000	40 08496 977178	800.0	1200.0	1190.0	160.200

Layer / Pallet	4	Master carton / Layer	4	Master carton / Pallet	16
Minimum pieces / shipment unit	10				

Generated on 09.01.2020, 10:29

All mentioned products and specifications are subject to change.
 In case of non-availability of any product or item, VARTA Consumer Batteries GmbH & Co. KGaA cannot be held responsible.
 Some of the products shown may not be available in all countries.
 The pictures might be different to the original products. Technical changes and errors reserved.

VARTA Consumer Batteries GmbH & Co. KGaA | Alfred-Krupp-Str. 9 | 73479 Ellwangen | Germany
 P.O. Box 1363 | 73473 Ellwangen | Germany | Phone (+49) 7961 83 0 | Fax (+49) 7961 83 440

